

2-494-**

2-495-**

2-496-**

基準分銅型円筒分銅寸法図

形状： 基準分銅型円筒分銅

材質： 非磁性ステンレス

単位：mm

表す量	D	H	h
20kg	142	213.5	146.5
10kg	113	169.5	116.5
5kg	90	136.8	91.8
2kg	66	99.4	68.4
1kg	53	76.0	54.0
500g	41	63.5	44.5
200g	31	45.7	31.7
100g	25	35.5	24.5
50g	20	28.3	19.3
20g	15	20.0	13.8
10g	13	14.7	8.7
5g	10	13.0	7.0
2g	8	8.4	4.4
1g	6	6.9	3.9

